

Applied Scope Warmer



A large teal geometric graphic consisting of a dark teal triangle on the left and a lighter teal triangle on the right, meeting at a diagonal line.

A SIMPLE SOLUTION FOR BETTER VISIBILITY

Simple to use
Minimises lens fogging

The Answer to Laparoscopic Fogging

Eliminates the need for chemical gels or wipes, retains heat throughout lengthy procedures and minimises the need for repeated scope cleansing and defogging.

Scope Warmer Seal

Retains heat throughout lengthy procedures

Scope Warmer Bottle

Minimises fogging by pre-warming the laparoscope to body temperature

Scope Warmer Base

Stabilises the system when a scope is inserted



PSI \leq 35
mmHg \leq 1810

APPLIED MEDICAL
Marino Santa Margherita, CA 90668
U.S.A.



Applied Scope Warmer

Reorder No.	Description	Qty/Box
C3001	Scope Warmer Bottle* Reusable Volume Capacity = 500mL	1 EACH
C3002	Scope Warmer Base* Reusable	1 EACH
C3101	Scope Warmer Seals Disposable, Single Use Individually packaged in sterile pouches	20/BOX

*Recommended Sterilization Method: Autoclave or STERRAD**100S system.
Please see Instructions for Use (IFU) for complete information.
Scope is not included.

Visit www.appliedmedical.com/laparoscopic
for more information

© 2018 Applied Medical Resources Corporation. All rights reserved.
Applied Medical and the Applied Medical logo design are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. **STERRAD is a U.S. registered trademark of Advanced Sterilization Products. SC01094B-EN

